

## 1. What is hexed?

The Hexed is short for Hexapod Daemon. The IRTF secondary mirror is mounted on a Physik Instrumente HEXEAPOD unit. For details about the HEXAPOD unit refer to the vendor's documentation. The HEXAPOD Controller (HEXC) provides serial data communication to the HEXAPOD unit.

The HEXC's RS-232 interface is connected to digimim, a Digi portserver II terminal server. digimim provides access to the HEXC RS-232 interface over TCP/IP. The Digi Portserver only services a single TCP/IP connect per port. The hexed daemon runs on an IRTF computer. Its purpose is to multiplex this single TCP/IP port, providing users and telescope application simultaneous access to the HEXC serial interface. Hexed replaces the Tip-Tilt ttconx/hexe programs that ran on mother (ttconx/hexe only provided exclusive access to HEXC's RS-232).

Hexed actually provide up to 5 concurrent TCP/IP connections. Normally host programs should connect, communication, then disconnect to the Hexed. Used in some cases permanent socket connection are used (Telnet for example).

## 2. Starting and Stopping the HexeD

Hexed is automatically started when the hexed\_host is running. To stop and restart the daemon you should used the start/stop script /etc/init.d/hexed, for example do:

1. Login to hexed\_host and become root.
2. Kill the hexed (if it is running):  

```
% /etc/init.d/hexed stop
```
3. Start the hexed  

```
% /etc/init.d/hexed start
```

### 3. How do I use HexeD?

Here are 3 methods to use the HexeD to command the HEXAPOD. They are presented here in order of preference:

#### 3.1 HexeGUI

An X based GUI to the HexeD is available on all IRTF Suns. HexeGUI is the HEXAPOD interface for observers and the IRTF staff. Refer to the HexeGUI Documentation.

#### 3.2 HexeIO

A simple hexeio program does a single command/reply sequence. This is a useful utility for 3<sup>rd</sup> party application and scripts. Refer to the HexeIO documentation.

#### 3.3. Telnet

You may also telnet to the hexed. Up to 5 network connections are provided by the hexed. A telnet session consumes 1 network connection. The user should exit when the telnet session is not needed. To open a telnet session you must be logged into one of the IRTF workstations on the summit (a firewall restriction). Use Stefan or boltzmann. Next, telnet to the hexed\_host using port 30015. There is an example session. The bold letters indicate what the user typed:

```
> telnet duke 30015
Trying 128.171.165.3...
Connected to duke.
Escape character is '^]'.
version
HexeDaemon v02.08 Aug 29 2002
Right now is 08/29/2002 12:48:56
no error -HD.EOC
stat n20
Commanded hexapod positions
Dimension          X          Y          Z          U          V          W
Transl. [mm]      -1.2500   -1.5000   -1.7010
Rotation [rad]                                0.0125    0.0067    0.0000
no error -HD.EOC
exit
closing connection.
Connection closed by foreign host.
```

After connecting with telnet, you can type HEXAPOD, or built in hexed commands.

Type 'exit' to close the session.

## 4. HexeD built-in commands

Normally the hexed just passes commands to the hexapod. There are a few internal hexed commands. These commands are described here:

**die** – Kills hexed.

Syntax `die`

**exit** – closes the telnet session used by the current client.

Syntax `exit`

**focus** – Set the focus value by setting the Z value on the hexapod.

Syntax `Focus z`

Range `z` – Value for Z. Limit to -8.0 to 8.0.

**getfocus** – prints the current Z values. Z indicated the focus value..

Syntax `GetFocus`

**help** – prints a summary of all the internal hexed commands.

Syntax `help`

**hexe** – passes a string to the hexed.

Syntax `Hexe string`

Range `string` – Any legal Hexapod Command.

**IncFocus** – Changes the focus value by incrementing the Z value on the hexapod by N.

Syntax `IncFocus n`

Range `n` – A value between -2.0 to 2.0. Note the resulting Z value is limited to -8.0 to 8.0.

**log** – Make the hexed writes a string to stdout.

Syntax `log string`

Range `string` – Any string.

**show.GV** – Prints out some hexed application variables.

Syntax `Show.GV`

**verbose** – Set the level of debug output printed by the daemon.

Syntax `verbose v`

Range `0 to 2` – The Lowest is 0 (default). Higher the value the more debug output.

**version** – Prints out the name and version number for hexed.

Syntax `version`

**sim** – Set the simulation flag. In simulation mode, no data is set to the digiport. Any HEXAPOD string are just reversed and returned to the caller.

Syntax `sim { off | on }`

Range `off` – disable simulation (default).

`on` – enables simulation..

**Z** – Set the focus value by setting the Z value on the hexapod. (same as Focus).

Syntax `z z`

Range `z` – Value for Z. Limit to -4.0 to 4.0.

## 5. HexeD FAQ

### 5.1 Where is the source code?

Contact the IRTF.

### 5.2 Any details on the RS-232 to digiPortServer II interface?

The HEXC's RS-232 was connected to port 2016 of a Digi PortServer II named digimim. The configuration for port #16 is

```
set line range=16 parity=N csize=8 error=null
set line range=16 baud=9600 stopb=1 break=ignore inpck=off istrip=off onlcr=off otab=off
set flow range=16 ixon=off aixon=off ixoff=off ixany=off itoss=off altpin=off
set flow range=16 rts=on dtr=off cts=on dcd=off dsr=off ri=off
set keys range=16 xon=^Q xoff=^S xona=^Q xoffa=^S
set keys range=16 eof=^D erase=^H intr=^C kill=^U tesc=^]
set port range=16 dev=prn sess=4 termtype=host uid=0 edelay=1 auto=off bin=off group=0
set port range=16 dest=0.0.0.0 dport=0
```

### 5.4 Prior the digimin. Hexapod was connected to digicass, a port server 16. What was the configuration information for that setup?

The HEXC's RS-232 was connected to port 2004 of a Digi PortServer 16 named digicass. The configuration for port #4 is:

```
set ports range=4 termtype="host" dev=prn sess=4 uid=0 edelay=1^M
set ports range=4 auto=off bin=off group=0 dport=0 dest=255.255.255.255^M
set line range=4 baud=9600 csize=8 parity=N stopb=1 break=ignore^M
set line range=4 error=null inpck=off istrip=off onlcr=off otab=off^M
set flow range=4 ixon=off aixon=off ixoff=off ixany=off itoss=off altpin=off^M
set flow range=4 rts=on dtr=off cts=on dcd=off dsr=off ri=off^M
set keys range=4 eof=^D erase=^H flush=^@ intr=^C kill=^E lnxt=^V^M
set keys range=4 tesc=^] xon=^S xoff=^Q xona=^S xoffa=^Q^M
set logins range=4 cmdprompt="digi>" logprompt="login:"^M
set logins range=4 passprompt="passwd:"^M
set logins range=4 write=off login=on passwd=on verbose=off^M
```

On Nov 14, 2002 the serial port was moved from digicass to digimim (a portserver II).

### 5.5 If I have a question about hexed, who should I contact?

The hexed was written by Tony Denault [denault@irtf.ifa.hawaii.edu](mailto:denault@irtf.ifa.hawaii.edu).